

1 **Bacterial starter cultures induce suitable changes in milk fatty acid profiles at different**
2 **fermentation conditions.**

3

4 Gholamreza Bahrami, Ali Mostafaie, Amir Kiani and Maryam Chalabi

5

6 SUPPLEMENTARY FILE

7

8

Table S1. The gas chromatograph conditions and the oven temperature.

Gas chromatograph conditions			The oven temperature		
Injector	Tem.	270 °C	Temp. (°C)	Rate (°C/min)	Hold (min)
	Split Ratio	1:70	50	-	1
Detector	Type	FID	100	10	2
	Tem.	300 °C	175	5	3
Column	Carrier Gas	Nitrogen	190	5	3
	Flow	1.5 mL/min	195	0.1	5
	CP-SIL 88,100 m×0.25mm×0.2µm		240	10	5

9

10